

Press Release

Facelift for reusable pipette tip boxes

Hamburg, August 2021

The new box is now called Box 2.0. The new design is a consistent, modern further development of the original product outfit.

When epT.I.P.S.[®] pipette tips were introduced by Eppendorf in 2002, they were already marketed in the popular reusable box. After 19 years, the epT.I.P.S. brand is now getting a new look.

But since looks aren't everything, the now even safer stackability of several boxes on top of each other is also ensured.

The new design also comes with changes suggested by customers. The viewing slits on the back of the box have largely disappeared, allaying customer concerns about possible contamination hazards.

Of course, the new Box 2.0 is also autoclavable up to 100 times and available in 3 sizes, suitable for epT.I.P.S. pipette tip sizes up to 5 mL - equivalent to 120 mm tip length.

More information on www.eppendorf.com/epTIPS-News

Press Contact:

Bettina Grochulla
Barkhausenweg 1
22339 Hamburg, Germany

Phone: +49 40 538 01 - 640
grochulla.b@eppendorf.de

www.eppendorf.com

About Eppendorf:

Eppendorf is a leading life science company that develops and sells instruments, consumables, and services for liquid handling, sample handling, and cell handling in laboratories worldwide. Its product range includes pipettes and automated pipetting systems, dispensers, centrifuges, mixers, spectrometers, and DNA amplification equipment as well as ultra-low temperature freezers, fermentors, bioreactors, CO2 incubators, shakers, and cell manipulation systems. Consumables such as pipette tips, test tubes, microplates, and single-use bioreactor vessels complement the range of highest-quality premium products.

Eppendorf products are most broadly used in academic and commercial research laboratories, e.g., in companies from the pharmaceutical and biotechnological as well as the chemical and food industries. They are also aimed at clinical and environmental analysis laboratories, forensics, and at industrial laboratories performing process analysis, production, and quality assurance.

Eppendorf was founded in Hamburg, Germany in 1945 and has about 4.500 employees worldwide. The company has subsidiaries in 26 countries and is represented in all other markets by distributors.

Press Contact:

Bettina Grochulla
Barkhausenweg 1
22339 Hamburg, Germany

Phone: +49 40 538 01 - 640
grochulla.b@eppendorf.de

www.eppendorf.com