

Press Release

Take Full Advantage of Versatility: Eppendorf SmartExtender now available

Hamburg, September 2018

The Eppendorf ThermoMixer® series is an established global standard for heating, cooling, and mixing of samples in the lab. All common vessels can be handled by specialized SmartBlocks on top of the mixer.

The new Eppendorf SmartExtender offers a comfortable incubation tool which can easily be used as an add-on to existing mixers and related SmartBlocks in the lab.

Up to 12x 1.5 mL vessels can be incubated in parallel whereas the temperature control for heating is independent from the SmartBlock in-use. Especially when performing a multi-step workflow like enzymatic reactions or needing a standard incubation temperature of e.g. 37°C in parallel to individual workflows, a second temperature within the same instrument saves time.

By taking advantage of the ThermoTop power support connectors, every existing Eppendorf ThermoMixer C, Fx, or ThermoStat C can be used with the SmartExtender.

For more information, please visit: www.eppendorf.com/thermomixer

Press Contact:

Bettina Doerler
Barkhausenweg 1
22339 Hamburg, Germany

Phone: +49 40 538 01-640
Fax: +49 40 538 01-662
doerler.b@eppendorf.de

www.eppendorf.com

About Eppendorf:

Eppendorf is a leading life science company that develops and sells instruments, consumables, and services for liquid handling, sample handling, and cell handling in laboratories worldwide. Its product range includes pipettes and automated pipetting systems, dispensers, centrifuges, mixers, spectrometers, and DNA amplification equipment as well as ultra-low temperature freezers, fermentors, bioreactors, CO2 incubators, shakers, and cell manipulation systems. Consumables such as pipette tips, test tubes, microplates, and single-use bioreactor vessels complement the range of highest-quality premium products.

Eppendorf products are most broadly used in academic and commercial research laboratories, e.g., in companies from the pharmaceutical and biotechnological as well as the chemical and food industries. They are also aimed at clinical and environmental analysis laboratories, forensics, and at industrial laboratories performing process analysis, production, and quality assurance.

Eppendorf was founded in Hamburg, Germany in 1945 and has more than 3,100 employees worldwide. The company has subsidiaries in 26 countries and is represented in all other markets by distributors.

Press Contact:

Bettina Doerler
Barkhausenweg 1
22339 Hamburg, Germany

Phone: +49 40 538 01-640
Fax: +49 40 538 01-662
doerler.b@eppendorf.de

www.eppendorf.com