

“Top Performance”

Varispenser 2 and Varispenser 2x –The new Uncompromising Solution for the most Demanding dispensing Applications

Hamburg, august 2017

Bottletop dispensers are needed in any lab where aggressive solutions like lyes, acids, bases or solvents are dispensed from large supply bottles. Bottletop dispensers need to be universally applicable and of uncompromising quality.

As a specialist for Liquid Handling laboratory products, Eppendorf has enhanced the proven quality and technology of the previous models in many details to develop the new models Varispenser 2 and Varispenser 2x. Thanks to the proven positive displacement principle, loss of residue is almost completely avoided. This technology combined with the high chemical resistance of all components that might come into contact with aggressive liquids allows the use of Varispenser 2 and Varispenser 2x with almost all liquids commonly used in labs. Accessories such as an easily assembled drying tube, a filter or a flexible discharge tube enable a broad range of dispensing options. All devices are equipped with standard GL 45 threads and can be used with almost all kinds of threads used in labs thanks to the additional adapters included in the scope of delivery.

The oval-shaped housing and the new design of the slider and the valve toggle of Varispenser 2x emphasize the modernity of the devices while ensuring optimal handling.

Press Contact:

Bettina Doerler
Barkhausenweg 1
22339 Hamburg

Tel: +49 40 538 01-640
Fax: +49 40 538 01-662
doerler.b@eppendorf.de

www.eppendorf.com

About Eppendorf AG:

Eppendorf is a leading life science company that develops and sells instruments, consumables, and services for liquid handling, sample handling, and cell handling in laboratories worldwide. Its product range includes pipettes and automated pipetting systems, dispensers, centrifuges, mixers, spectrometers, and DNA amplification equipment as well as ultra-low temperature freezers, fermentors, bioreactors, CO2 incubators, shakers, and cell manipulation systems. Consumables such as pipette tips, test tubes, microliter plates, and single-use bioreactor vessels complement the range of highest-quality premium products.

Eppendorf products are most broadly used in academic and commercial research laboratories, e.g., in companies from the pharmaceutical and biotechnological as well as the chemical and food industries. They are also aimed at clinical and environmental analysis laboratories, forensics, and at industrial laboratories performing process analysis, production, and quality assurance.

Eppendorf was founded in Hamburg, Germany in 1945 and has more than 3,000 employees worldwide. The company has subsidiaries in 26 countries and is represented in all other markets by distributors.

Press Contact:

Bettina Doerler
Barkhausenweg 1
22339 Hamburg

Tel: +49 40 538 01-640
Fax: +49 40 538 01-662
doerler.b@eppendorf.de

www.eppendorf.com